

GERSTEL AppNote 159

Thermal Gravimetric Analysis/Mass Spectrometry Simulation using the GERSTEL PYRO Core System

John R. Stuff¹, Jackie A. Whitecavage¹, Michael McAdams², Yunyun Nie³, and Eike Kleine-Benne³

¹Gerstel, Inc., 701 Digital Dr. Suite J, Linthicum, MD 21090, USA

²Adhesives Research, 400 Seaks Run Road, Glen Rock, PA, 17327, USA

³Gerstel GmbH & Co. KG, Eberhard-Gerstel-Platz 1, 45473 Mülheim an der Ruhr, Germany

Keywords

Pyrolysis, thermogravimetric analysis, gas chromatography, mass spectrometry, evolved gas analysis

Abstract

This study shows the use of the GERSTEL PYRO Core system for generating simulated thermal gravimetric-mass spectrometry data (TGA-MS) of polymer samples. This mode of operation is also referred to as evolved gas analysis (EGA).

For EGA analysis, a short piece of uncoated capillary is attached from the GC inlet to the mass spectrometer. A relatively slow temperature ramp is run similar to that used for actual TGA analysis. The total ion chromatogram can be plotted as a function of temperature to produce simulated TGA-MS data. The mass spectral data is used to identify degradation products at various temperatures.

Several types of polymers are examined in this study. Simulated TGA-MS data is compared with actual TGA data.

Introduction

The GERSTEL PYRO Core system with its integrated MultiPurpose Sampler (MPS), pyrolyzer, and Cooled Injection System (CIS 4), PTV-type inlet, provides the user with a multitude of analytical options to utilize for sample analysis. It offers efficient automation (up to 2 x 98 samples) and a variety of pyrolysis modes including standard pulsed, sequential, and fractionated pyrolysis. The CIS 4 can be used to cryofocus analytes in the inlet or as a hot split interface for direct transfer to the column.

This study describes the use of the GERSTEL PYRO Core system for simulation of thermogravimetric (TGA)/MS analysis. The EGA mode relies on a relatively slow temperature ramp for the pyrolyzer and a short piece of deactivated fused silica to connect the CIS to the mass spectrometer.

The thermogravimetric data was collected on a TA Instruments Model Q5000 thermogravimetric analyzer.

Several polymer types were analyzed in order to demonstrate this mode with the GERSTEL PYRO Core system.

GERSTEL AppNote 159

Experimental

Instrumentation

Analyses were performed using a GERSTEL PYRO Core system on a 7890 equipped with a 5975C Inert XL MSD with triple axis detector (Agilent Technologies).

Sample Preparation - Pyrolysis

Approximately one to three milligrams of sample were weighed into a short, quartz test-tube shaped pyrolysis vessel. The quartz tubes were connected to pyrolysis adapters and placed into a 98 position pyrolysis tray in the MPS.

Sample Preparation - TGA Analysis

Approximately 20 mg of sample was placed in the TGA pan.

Analysis Conditions PYRO Core System

Pyrolysis Lead Time 0.10 min
 Follow up Time 0.50 min
 Initial Time 1.5 min
 80 °C (1.0 min); 0.17 °C/s to 600 °C (1.0 min)
 solvent vent (50 mL/min)
 -120 °C; 12 °C/min to 280 °C (3 min)

CIS 4

Analysis Conditions Agilent 7890 GC

Column 1 m Deactivated Fused Silica (Restek)
 $d_i = 0.1 \text{ mm}$
 Pneumatics He, constant flow = 1 mL/min
 $p_i = 7.1 \text{ psi}$
 Oven 250 °C (isothermal)

Analysis Conditions TA Instruments Q5000

TGA 50 °C; 20 °C/min to 700 °C;
 100 °C/min to 1000 °C

Results and Discussion

Four polymer samples were used to demonstrate the GERSTEL PYRO Core system used in EGA mode. Two Kraton® polymers (styrenic block copolymers consisting of polystyrene and rubber blocks), a urethane elastomer, and an acrylonitrile-styrene-acrylate co-polymer.

Figure 1 shows the TGA and first derivative curves for a Kraton® D1161 sample. The curve shows rapid weight loss starting at 325 °C and continuing up to 400 °C. A second transition is seen in the first derivative curve from 415 °C to 470 °C. Spikes seen in the TGA curves from 370-400 °C, are caused by balance being perturbed by the rapid weight loss of the sample during this time frame.




Figure 1: TGA and first derivative curve for Kraton® D1161.

GERSTEL AppNote 159

Figure 2 shows the TGA-MS data collected using the GERSTEL pyrolyzer in EGA mode. The shape is similar to that seen in the first derivative TGA curve. The spectral information obtained shows

that the major, first, transition is mainly from loss of the natural rubber component of the polymer. The second smaller transition, 415-500 °C, is mainly from loss of the polystyrene.




Figure 2: TGA-MS data for Kraton® D1161

Figure 3 shows the TGA and first derivative curves for a Kraton® D1111 sample. The curve shows rapid weight loss starting at 325 °C and continuing up to 387 °C. A second transition is seen in the first derivative curve at 434 °C. In comparison to the other

Kraton® sample (Figure 1), the main transition peaks at a lower temperature and the second transition is a larger percentage of the overall weight.




Figure 3: TGA and first derivative curve for Kraton® D1111.

GERSTEL AppNote 159

The TGA-MS data, Figure 4, shows that the decomposition products for the two transitions are similar. The spectral information obtained shows that the major, first, transition is mainly from loss of

the natural rubber component of the polymer. The second smaller transition is mainly from loss of the polystyrene.




Figure 4: TGA and first derivative curve for Kraton® D1111.

Figure 5 shows the TGA and first derivative curves for a urethane elastomer. The first derivative curve shows two maxima at 389 °C

and 432 °C. The TGA-MS total ion chromatogram (not shown) shows the same two peaks.




Figure 5: TGA and first derivative curve for a urethane elastomer.

GERSTEL AppNote 159

Figure 6 shows Extracted Ion Chromatograms (EICs) for the urethane elastomer. The EICs show that the initial loss is due to loss of the monomer, diphenylmethane diisocyanate ((MDI), $m/z=250$)

and the linking agent butanediol ($m/z=42$). The second maximum is caused by loss of MDI-butanediol ($m/z=129$).




Figure 6: TGA-MS EICs for a urethane elastomer.

Figure 7 shows the TGA and first derivative curves for an acrylonitrile-styrene-acrylate copolymer. The first derivative curve shows a

single maximum at 442 °C. The TGA-MS total ion chromatogram (not shown) also shows a single peak.




Figure 7: TGA and first derivative curve for an acrylonitrile-styrene-acrylate copolymer.

GERSTEL AppNote 159

When specific ion chromatograms are extracted from the TIC, as shown in figure 8, a clearer picture of the polymer degradation develops. The EICs show that around 325 °C, the ions for styrene ($m/z=104$) and butanol ($m/z=56$) appear. The butanol loss reaches

a maximum at 400 °C, and ends around 425 °C. The styrene loss reaches a maximum around 430 °C along with the ion for the acrylonitrile monomer ($m/z=53$) and the ion for butylmethacrylate ($m/z=73$).




Figure 8: TGA-MS EICs for an acrylonitrile-styrene-acrylate copolymer.

Conclusions

The GERSTEL PYRO Core system can be used for evolved gas analysis, producing data similar to thermal gravimetric analysis-mass spectrometry. This type of data in conjunction with standard TGA analysis can provide useful information regarding the thermal stability and degradation pathway for polymers as well as subtle differences between polymers of similar composition. This type of data can also be used to determine appropriate polymer pyrolysis temperatures.