

GERSTEL AppNote 202

Analysis of Food Samples using Thin Film Solid Phase Microextraction (TF-SPME) and Thermal Desorption GC/MS

Laurel Vernarelli, Jackie Whitecavage, John Stuff

GERSTEL, Inc., 701 Digital Drive, Suite J, Linthicum, MD 21090, USA

Keywords

LabWorks Platform, Thin Film Solid Phase Micro-Extraction (TF-SPME), Gas Chromatography, Mass Spectrometry, Flavor Compounds

Abstract

Food products are routinely monitored for quality, authenticity, and safety. Analysis of food products may also be necessary due to off flavor or odor complaints. The aroma and flavor profiles of each product are unique and made up from a variety of semi-volatile and volatile compounds including aldehydes, ketones, acids, alcohols, terpenes, esters and other trace level components. The wide range of concentrations, polarities, and functional groups used in the composition of a flavor/aroma profile can make the analysis of the sample difficult. Techniques which are simple, use little or no solvent and encompass a wide range of analytes are desirable.

TF-SPME is an extension of regular SPME. It is more sensitive than regular SPME due to the increased surface area and phase volume, both of which lead to improved analyte recovery. The TF-SPME device is a 20 mm x 4.8 mm carbon mesh sheet impreg-

nated with a sorptive phase. The TF-SPME devices can be used in headspace or immersive mode.

In headspace mode, the TF-SPME device is suspended above a solid or liquid sample inside a sealed sample vial. In immersive mode, it is placed directly in a liquid sample. In both cases, the sample is agitated by stirring. After extraction, the devices are blotted dry and placed in an empty thermal desorption tube. They are then analyzed by thermal desorption GC/MS, a technique that is easy to automate, simple to use and requires no solvent.

The work presented here demonstrates the application of Thin Film Solid Phase Microextraction (TF-SPME) to the determination of aroma and flavor components in a variety of high fat food products.

GERSTEL AppNote 202

Introduction

The GERSTEL LabWorks Platform provides the user with a multitude of analytical options to utilize for chemical analysis. This study describes the use of the GERSTEL MPS in combination with a CIS 4 inlet and TDU for automated thermal desorption of TF-SPME devices with a divinylbenzene/polydimethylsiloxane coating (DVB/PDMS) that have been used to extract flavor components from a variety of food samples. Among these are: blue cheese, dark chocolate, Caesar dressing, and strawberry cream cheese. Monitoring of key flavor components is important for quality control and product development of foods. Identification of off-aromas can help pinpoint and find a remedy for product defects.

Experimental

Instrumentation

GERSTEL LabWorks Platform with Agitator/Stirrer on Agilent 7890A GC / 5977B MSD

Analysis Conditions LabWorks Platform

TDU splitless
40 °C (0.5 min); 720 °C/min to 280 °C (3 min)

PTV solvent vent (50 mL/min), split 5:1
-150 °C; 12 °C/s to 280 °C (3 min)

Analysis Conditions GC

Agilent 7890A GC

Dark chocolate:

Pneumatics He, constant flow, 1 mL/min
Column 30 m Rxi-5Sil MS (Restek)
 $d_i = 0.25 \text{ mm}$ $d_f = 0.25 \text{ }\mu\text{m}$
Oven 40 °C (3 min), 7 °C/min to 280 °C (3 min)

Other samples:

Pneumatics He, constant pressure, 336 kPa
Column 30 m DB-Wax (Agilent)
 $d_i = 0.25 \text{ mm}$ $d_f = 0.25 \text{ }\mu\text{m}$
Oven: 40 °C (1 min), 10 °C/min to 230 °C (5 min)

Sample Preparation

All analyzed food products were purchased at a local store. A 2.0 g aliquot of each sample was weighed and placed in an empty 20 mL screw cap vial. A TF-SPME device was suspended in the head-space of the vial. The vial was placed into the GERSTEL agitator/stirrer and extracted at 40 °C for 3 hours at 250-500 rpm agitation speed. The strawberry cream cheese sample was also extracted at room temperature (25 °C) for 15 hours.

Results and Discussion

All food samples, apart from dark chocolate, were separated using a DB WAX column. The dark chocolate food sample was separated using a Rxi-5Sil MS column.

GERSTEL AppNote 202

Dark Chocolate

TF-SPME was applied to the headspace of the dark chocolate sample. Figure 1 shows the resulting total ion chromatogram with identified peaks labeled. A wide variety of volatile and semi volatile flavor compounds were observed including pyrazines, esters, alcohols, acids, aldehydes, ketones, furans, etc. Several pyrazine compounds, which impart important cocoa/nutty flavors to the chocolate, were observed: 2,5-dimethyl-2-(3-methylbutyl)-pyrazine, 3-ethyl-2,5-dimethyl-pyrazine, tetramethyl pyrazine (roasty

note) and 2,3,5-trimethyl-6-ethylpyrazine. Other key aroma compounds present in large quantities include acetic acid, isovaleraldehyde, 2-methylbutyraldehyde, 2,3-butandiol, tetrahydrofurfuryl alcohol and alpha-pinene.

The results show that TF-SPME is a highly useful technique for assessing off-aromas or quality defects in dark chocolate with very little sample preparation.




Figure 1: Total ion chromatogram resulting from Thin Film Solid Phase Microextraction of dark chocolate.

GERSTEL AppNote 202

Blue Cheese

TF-SPME was applied to the headspace of a blue cheese sample. Figure 2 shows the resulting total ion chromatogram with identified peaks labeled. A wide variety of compound types were observed including aldehydes, alcohols, esters, acids and phenols. The predominant compounds are 2-pentanone, 2-heptanone and 2-nonanone, which are metabolites formed by *Penicillium* mold. These methyl ketones are responsible for the characteristic fla-

vor of mold-ripened cheese. Other flavor compounds that were observed in the chromatogram include the secondary alcohols 2-pentanol, 2-heptanol, 2-nonanol, as well as the fatty acids butanoic acid, hexanoic acid, octanoic acid, nonanoic acid and n-decanoic acid. These results show that TF-SPME is a highly useful technique for assessing aroma compounds in blue cheese.




Figure 2: Total ion chromatogram resulting from Thin Film Solid Phase Microextraction of blue cheese.

GERSTEL AppNote 202

Caesar Dressing

TF-SPME was applied to the headspace of a caesar dressing sample. Figure 3 shows the resulting total ion chromatogram with identified peaks labeled. A wide variety of compounds were observed including terpenes/terpenoids, alcohols, acids and aldehydes. The predominant compounds are sorbic acid, which is a food preservative and ethyl acetate, a solvent used for flavors. A wide variety of terpenes and terpenoids are also found in the chromatogram, due to the presence of black pepper in the dressing.

These include pinene, limonene, phellandrene, and caryophyllene, key aroma compounds found in the outer fruit layer of black pepper. Fatty acids such as butanoic acid, hexanoic acid, heptanoic acid, octanoic acid and n-decanoic acid are also present, likely due to the presence of parmesan and romano cheese listed as ingredients in the dressing. These results show that TF-SPME is a highly useful technique for assessing off-aroma and quality defects in dressing samples.




Figure 3: Total ion chromatogram resulting from Thin Film Solid Phase Microextraction of caesar dressing.

GERSTEL AppNote 202

Strawberry Cream Cheese

TF-SPME was applied to the headspace of a strawberry flavored cream cheese sample. Figure 4 shows a chromatogram resulting from an extraction at 25 °C that lasted 16 hours. A wide variety of compounds were observed including alcohols, esters, aldehydes, acids, furans, and lactones. The largest peaks in the chromatogram

were identified as ethyl acetate, a flavor solvent; ethyl butyrate, a commonly used artificial flavoring for fruit-flavored foods; ethyl 2-methyl butyrate, a fruity flavoring agent; and acetoin, a food flavoring agent with a pleasant, buttery aroma.




Figure 4: Total ion chromatogram resulting from Thin Film Solid Phase Microextraction of strawberry flavored cream cheese, extraction at 25 °C for 15 hours.

GERSTEL AppNote 202

Figure 5 shows the resulting chromatogram following extraction at 40 °C for 3 hours. This chromatogram showed an increase in peak areas corresponding to Maillard reaction products: furfural, 2-furanmethanol, maltol, 3,5-dihydroxy-6-methyl-2,3-dihydro-4H-pyran-4-one and 5-hydroxymethylfurfural. This indicates that TF-SPME is well suited for monitoring undesirable Maillard reaction products that may occur due to food processing or heating of

food products, such as strawberry cream cheese.

The results show that TF-SPME is well suited for quantification as well as assessment of quality defects that may arise due to food processing or packaging, while requiring very little sample preparation.




Figure 5: Total ion chromatogram resulting from Thin Film Solid Phase Microextraction of strawberry flavored cream cheese, extraction at 40 °C for 3 hours.

Conclusions

The GERSTEL LabWorks Platform is a versatile platform for qualitative analysis of a wide variety of food products. Thin Film Solid Phase Microextraction (TF-SPME) is a robust technique that allows the user to quickly identify and quantify a wide variety of compound types in food samples. In addition to TF-SPME, the versatile GERSTEL MPS/TDU/CIS system allows several techniques to be applied quickly to the same samples, rapidly providing the necessary data for quality control, product development, troubleshooting, or competitive analysis. With the appropriate use

of standards, any of the techniques can be used for quantitative analysis.