

The GERSTEL logo consists of the word "GERSTEL" in a bold, black, sans-serif font, enclosed within a red rectangular border.

MAKING LABS WORK

TDU to TDU 2 upgrade  
until December 31, 2023

A detailed view of a GERSTEL TDU 2 instrument. It features a stainless steel upper section with a glass tube and a red and white lower section. The brand name "GERSTEL" is printed on the red part, and "TDU 2" is printed in large white letters on the red background. The instrument is mounted on a white surface.

GERSTEL

TDU  
2

A red circular badge with white text. The text reads "SPECIAL OFFER" in a smaller font, followed by "50%" in a very large font, and "off" in a smaller font below it.

SPECIAL  
OFFER  
50%  
off

The GERSTEL Thermal Desorption Unit TDU, a world-wide trendsetter and runaway success, was discontinued In 2017 and replaced by the TDU 2. The seven year clock until end of support runs until early 2024. At that time, spare parts will no longer be available and product support will end.

Upgrading to the TDU 2 is therefore urgent – and it brings great benefits!

TDU 2 is a 4<sup>th</sup> generation GERSTEL Thermal Desorber based on more than 20 years of experience.

The TDU 2 offers improved robustness and serviceability:

- Built-in alignment support for easy liner change
- New and improved transfer heater
- New and improved coolant and gas connections
- Improved, user friendly maintenance access
- Improved robustness, system monitoring, and serviceability

The TDU 2 enables ultra-trace analysis of a wide range of samples. Mounted directly on top of the CIS inlet without valve or transfer line in the sample flow path, the TDU 2 offers extreme inertness, best possible analyte recovery for accurate results, and lowest LODs.

**GERSTEL is offering a significant 50 % discount on the TDU to TDU 2 upgrade.**

**This special offer expires on December 31, 2023, so act now!**

To upgrade your TDU to a TDU 2, please contact your GERSTEL representative as soon as possible, or send a mail to: [info@gerstel.com](mailto:info@gerstel.com) with the subject line: **TDU Upgrade**

Please note: Service installation required.



**GERSTEL**

MAKING LABS WORK

GERSTEL GmbH & Co. KG  
Eberhard-Gerstel-Platz 1  
45473 Mülheim an der Ruhr  
Germany

[www.gerstel.com](http://www.gerstel.com)



GERSTEL®, GRAPHPACK®, TWISTER® and TWICESTER® are registered trademarks of GERSTEL GmbH & Co. KG. · Subject to change. Agilent® is a registered Trademark of Agilent Technologies, Inc. · Copyright by GERSTEL GmbH & Co. KG



 **Agilent Technologies**  
Premier Solution Partner